

# Anti-FCN1/FCNA antibody



Catalog Number: 106136

## Product name

Anti-FCN1/FCNA antibody

## Immunogen

[Human FCN1/FCNA \(His Tag\) recombinant protein](#)

## Specificity

Human FCN1

## Antibody description

Rabbit monoclonal to FCN1/FCNA

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human FCN1 (rh FCN1; NP\_001994.2; Met1-Ala326).

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions


IHC-P: 1-10 µg/mL

## Validations



Ficolin-1 / Ficolin-A / FCN1 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human FCN1 in human placenta with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to peripheral blood leukocytes.



Ficolin-1 / Ficolin-A / FCN1 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human FCN1 in human spleen (from 2 donors) with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to peripheral blood leukocytes.