

Anti-CD226/DNAM-1 antibody



Catalog Number: 106228

Product name

Anti-CD226/DNAM-1 antibody

Immunogen

[Mouse CD226/DNAM-1 \(His Tag\) recombinant protein](#)

Specificity

Mouse CD226 / DNAM-1

Antibody description

Rabbit monoclonal to CD226/DNAM-1

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD226 / DNAM-1 (rM CD226 / DNAM-1; NP_848802.2; Met1-Pro254) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be

disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



Ly108 / SLAMF6 Antibody, Rabbit MAb, Immunohistochemistry

CD226 / DNAM-1 Antibody (FITC) , Rabbit MAb, Flow cytometric analysis