

# Anti-CD226/DNAM-1 antibody



Catalog Number: 106228

## Product name

Anti-CD226/DNAM-1 antibody

## Immunogen

[Mouse CD226/DNAM-1 \(His Tag\) recombinant protein](#)

## Specificity

Mouse CD226 / DNAM-1

## Antibody description

Rabbit monoclonal to CD226/DNAM-1

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD226 / DNAM-1 (rM CD226 / DNAM-1; NP\_848802.2; Met1-Pro254) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be

disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

FCM

## Dilutions

FCM

## Validations



Ly108 / SLAMF6 Antibody, Rabbit MAb,  
Immunohistochemistry

CD226 / DNAM-1 Antibody (FITC) , Rabbit MAb,  
Flow cytometric analysis