

Anti-IL-1RA/IL1RN antibody



Catalog Number: 106362

Product name

Anti-IL-1RA/IL1RN antibody

Specificity

Human IL-1RA / IL1RN

Antibody description

Rabbit monoclonal to IL-1RA/IL1RN

Preparation

This antibody was obtained from a rabbit immunized with a synthetic peptide corresponding to the center region of the human IL1Ra / IL1RN.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 5-20 µg/mL

Validations



IL-1RA / IL1RN Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human IL1RA in human skin with rabbit monoclonal antibody (10 µg/mL, formalin-fixed paraffin embedded



sections).

IL-1RA / IL1RN Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human IL1RA in human kidney with rabbit monoclonal antibody (10 µg/mL, formalin-fixed paraffin embedded sections).