

Anti-Coagulation Factor X/F10 antibody



Catalog Number: 106398

Product name

Anti-Coagulation Factor X/F10 antibody

Immunogen

[Human Coagulation Factor X/F10 \(His Tag\) recombinant protein](#)

Specificity

Human Coagulation Factor X / F10

Antibody description

Mouse monoclonal to Coagulation Factor X/F10

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human Coagulation Factor X / F10 (rh Coagulation Factor X / F10; NP_000495.1; Met1-Lys488). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

IHC-P

Dilutions

IHC-P: 10-30 µg/mL

Validations



Coagulation Factor X / F10 Antibody, Mouse MAb, Immunohistochemistry

Immunochemical staining of human FX in human cirrhosis with mouse monoclonal antibody (15 µg/mL, formalin-fixed paraffin embedded sections).