

Anti-MYC antibody



Catalog Number: 106415

Product name

Anti-MYC antibody

Specificity

Recognize N-terminal, C-terminal in fusion proteins.

Antibody description

Mouse monoclonal to MYC

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, a synthetic Myc tag peptide conjugated to KLH. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

IP

Dilutions

IP: 1-4 μ g/mg of lysate

Validations



Items	Lanes	A	B
Sample (whole cell lysate)		HA-ARG1-myc	myc-ARG1-HA
Sample quantity			0.5 mg
IP antibody quantity		2 μ g [Cat# 100029-M007]	
Protein G agarose		15 μ l of 50% Protein G agarose	
Gel		13% SDS-PAGE reducing gel	
Primary antibody		Myc tag antibody at 1 μ g/ml [Cat# 100029-M007]	
Secondary antibody		DyLight 800-labeled antibody to mouse IgG (H+L), at 1:7500 dilution	

Anti-Myc Tag Antibody, Mouse MAb, Immunoprecipitation