

# Anti-Syndecan-4/SDC4 antibody



Catalog Number: 106436

## Product name

Anti-Syndecan-4/SDC4 antibody

## Immunogen

[Mouse Syndecan-4/SDC4 \(His Tag\) recombinant protein](#)

## Specificity

Mouse Syndecan-4 / SDC4

**No cross-reactivity** in ELISA with Mouse SDC1

## Antibody description

Rabbit monoclonal to Syndecan-4/SDC4

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse Syndecan-4 / SDC4 (rM Syndecan-4 / SDC4; O35988; Met1-Val146).

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions


IHC-P: 1-10 µg/mL

## Validations



Syndecan-4 / SDC4 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse SDC4 in mouse stomach with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded



sections).

Syndecan-4 / SDC4 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse SDC4 in mouse intestine with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).