

# Anti-TROP-2/TACSTD2 antibody



Catalog Number: 106528

## Product name

Anti-TROP-2/TACSTD2 antibody

## Immunogen

[Mouse TROP-2/TACSTD2 \(His Tag\) recombinant protein](#)

## Specificity

Mouse TROP-2 / TACSTD2

**Has cross-reactivity** in ELISA with Human TROP2

## Antibody description

Rabbit monoclonal to TROP-2/TACSTD2

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse TROP-2 / TACSTD2 (rM TROP-2 / TACSTD2; Q8BGV3; Met1-Gln270).

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions


IHC-P: 0.05-2 µg/mL

## Validations



TROP-2 / TACSTD2 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse TACSTD2 in mouse kidney with rabbit monoclonal antibody (0.05 µg/mL, formalin-fixed paraffin embedded



sections).

TROP-2 / TACSTD2 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse TACSTD2 in mouse skin with rabbit monoclonal antibody (0.05 µg/mL, formalin-fixed paraffin embedded sections).