

Anti-CD36/SCARB3 antibody



Catalog Number: 106681

Product name

Anti-CD36/SCARB3 antibody

Immunogen

[Rat CD36/SCARB3 \(His Tag\) recombinant protein](#)

Specificity

Rat CD36 / SCARB3

Antibody description

Rabbit monoclonal to CD36/SCARB3

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rat CD36 / SCARB3 (rR CD36 / SCARB3; Q66HI4; Gly30-Lys439).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 1-10 µg/mL

Validations



CD36 / SCARB3 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of rat CD36 in rat heart with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).



CD36 / SCARB3 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of rat CD36 in rat small intestine with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).