

# Anti-p53/TP53 antibody



Catalog Number: 106708

## Product name

Anti-p53/TP53 antibody

## Specificity

human P53 phosphorylated at Ser392

## Antibody description

Rabbit polyclonal to p53/TP53

## Preparation

Produced in rabbits immunized with A synthetic phospho-peptide corresponding to residues surrounding Ser392 of human P53., and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 10-30 µg/ml

## Validations



Lanes	A
Items	
Sample (whole cell lysate)	A431
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended concentration	1.5 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody To Mouse IgG (H+), at 1:7500 dilution.
Developed using Odyssey imaging system.	
Explanation	Predicted band size : 52 kDa Observed band size : 52 kDa

P53 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot