

Anti-Neurexophilin-1/NXPH1 antibody



Catalog Number: 106742

Product name

Anti-Neurexophilin-1/NXPH1 antibody

Immunogen

[Rat Neurexophilin-1/NXPH1 \(His Tag\) recombinant protein](#)

Specificity

Rat Neurexophilin-1 / NXPH1

Antibody description

Rabbit monoclonal to Neurexophilin-1/NXPH1

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rat Neurexophilin-1 / NXPH1 (rR Neurexophilin-1 / NXPH1; Q63366-1; Met1-Gly271).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 0.05-2 µg/mL

Validations



Neurexophilin-1 / NXPH1 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of rat NXPH1 in rat brain with rabbit monoclonal antibody (0.05 µg/mL, formalin-fixed paraffin embedded



sections).

Neurexophilin-1 / NXPH1 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of rat NXPH1 in rat spinal cord with rabbit monoclonal antibody (0.05 µg/mL, formalin-fixed paraffin embedded sections).