

PSEN2-Specific antibody

Catalog Number: 106963

Product name

PSEN2-Specific antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PSEN2-Specific Rabbit Polyclonal antibody.
Positive IHC detected in human kidney tissue,
human breast cancer tissue, human colon tissue,
human placenta tissue.

Preparation

This antibody was obtained by immunization of
Peptide (Accession Number: NM_000447).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded
human kidney tissue slide using Catalog
No:114300(PSEN2-Specific Antibody) at dilution
of 1:200 (under 10x lens)



Immunohistochemistry of paraffin-embedded
human kidney tissue slide using Catalog
No:114300(PSEN2-Specific Antibody) at dilution
of 1:200 (under 40x lens)