

FMO2-specific antibody



Catalog Number: 106972

Product name

FMO2-specific antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

FMO2-specific Rabbit Polyclonal antibody.
Positive WB detected in mouse kidney tissue,
mouse lung tissue. Observed molecular weight by
Western-blot: 50kd, 57kd

Preparation

This antibody was obtained by immunization of
Peptide (Accession Number: NM_001460).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse kidney tissue were subjected to SDS PAGE
followed by western blot with Catalog
No:110729(FMO2-specific antibody) at dilution of
1:300