

MEF2C-Specific antibody



Catalog Number: 106992

Product name

MEF2C-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MEF2C-Specific Rabbit Polyclonal antibody.
Positive IHC detected in human testis tissue, human brain tissue, human kidney tissue, human lung tissue, human spleen tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002397).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human testis using Catalog No:112637(MEF2C-Specific antibody) at dilution of 1:50 (under 10x lens)