

CXCR3B-specific antibody



Catalog Number: 107005

Product name

CXCR3B-specific antibody

Specificity

Human; other species not tested.

Antibody description

CXCR3B-specific Mouse Monoclonal antibody. Positive WB detected in HeLa cells, Jurkat cells, U-937 cells. Observed molecular weight by Western-blot: 40kd+43kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: AF469635). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107186(CXCR3B antibody) at dilution of 1:1000