

BID antibody

Catalog Number: 107018

Product name

BID antibody

Specificity

Human; other species not tested.

Antibody description

BID Mouse Monoclonal antibody. Positive WB detected in Jurkat cells, HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of BID recombinant protein (Accession Number: NM_001244567). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107018(BID



Antibody) at dilution of 1:1000

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107018(BID Antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells using Catalog No:107018(BID Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)