

BRCA1 antibody



Catalog Number: 107021

Product name

BRCA1 antibody

Specificity

Human; other species not tested.

Antibody description

BRCA1 Mouse Monoclonal antibody. Positive IHC detected in human breast cancer tissue, human ovary tumor tissue, human testis tissue. Positive WB detected in SKOV-3 cells. Observed molecular weight by Western-blot: 207 kDa

Preparation

This antibody was obtained by immunization of BRCA1 recombinant protein (Accession Number: BC115037). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



SKOV-3 cells were subjected to SDS PAGE followed by western blot with (BRCA1 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using (BRCA1 Antibody) at dilution of 1:100 (under 40x lens)