

CD138/Syndecan-1 antibody



Catalog Number: 107042

Product name

CD138/Syndecan-1 antibody

Specificity

Human; other species not tested.

Antibody description

CD138/Syndecan-1 Mouse Monoclonal antibody. Positive WB detected in BxPC-3 cells. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of CD138/Syndecan-1 recombinant protein (Accession Number: NM_002997). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107042(SDC1,CD138 antibody) at dilution of 1:400