

CD19 antibody



Catalog Number: 107046

Product name

CD19 antibody

Specificity

Human; other species not tested.

Antibody description

CD19 Mouse Monoclonal antibody. Positive WB detected in Raji cells, Daudi cells, Ramos cells. Observed molecular weight by Western-blot: 95 kDa

Preparation

This antibody was obtained by immunization of CD19 recombinant protein (Accession Number: NM_001770). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:107046(CD19 Antibody) at dilution of 1:2000