

# CD1A antibody



Catalog Number: 107047

---

## Product name

CD1A antibody

## Specificity

Human; other species not tested.

## Antibody description

CD1A Mouse Monoclonal antibody. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 49 kDa

## Preparation

This antibody was obtained by immunization of CD1A recombinant protein (Accession Number: BC031645). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:107047(CD1A antibody) at dilution of 1:300