

Antithrombin III antibody



Catalog Number: 107060

Product name

Antithrombin III antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Antithrombin III Mouse Monoclonal antibody.
Positive WB detected in L02 cells, human brain tissue, human testis tissue, rat brain tissue.
Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of Antithrombin III recombinant protein (Accession Number: BC022309). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:107060(SERPINC1 Antibody) at dilution of 1:1000