

APOD antibody



Catalog Number: 107064

Product name

Monoclonal

APOD antibody

Ig Type

Specificity

Mouse IgG1

Human; other species not tested.

Applications

Antibody description

ELISA, WB

APOD Mouse Monoclonal antibody. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 30-33 kDa

Dilutions

Recommended Dilution:


Preparation

WB: 1:5000-1:50000

This antibody was obtained by immunization of APOD recombinant protein (Accession Number: NM_001647). Purification method: Protein G purified.

Validations

Formulation



PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:107064(APOD Antibody) at dilution of 1:20000