

# 6\*His, His-Tag antibody

Catalog Number: 107093

## Product name

6\*His, His-Tag antibody

## Specificity

Tag

## Antibody description

6\*His, His-Tag Mouse Monoclonal antibody.  
Positive IP detected in Transfected HEK-293 cells.  
Positive WB detected in recombinant protein.  
Positive IF detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of Peptide. Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IP, IF

## Dilutions


Recommended Dilution:

WB: 1:5000-1:50000


IP: 1:5000-1:50000

IF: 1:10-1:100

## Validations




Western blot of 6\*His-tagged fusion protein with anti-6\*HIS tag (Catalog No:117294) at various dilutions.



dilutions.

IP Result of anti-6\*His, His-Tag (IP:Catalog No:117294, 7ug; Detection:Catalog No:117294 1:10000) with Transfected HEK-293 cells lysate



300ug.

Immunofluorescent analysis of Transfected HEK-293 cells using Catalog No:117294(6\*His, His-Tag Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)