

HRP-conjugated 6*His, His-Tag antibody



Catalog Number: 107094

Product name

Monoclonal

HRP-conjugated 6*His, His-Tag antibody

Ig Type

Specificity

Mouse IgG1

Tag

Applications

Antibody description

ELISA, WB

HRP-conjugated 6*His, His-Tag Mouse Monoclonal antibody. Positive WB detected in Recombinant protein. Observed molecular weight by Western-blot: -

Dilutions


Recommended Dilution:

WB: 1:5000-1:50000

Preparation

This antibody was obtained by immunization of 6xHis Peptide. Purification method: Protein G Purified.

Validations



Formulation

PBS with 0.02% Thimerosal, 50% Glycerol, pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Recombinant protein were subjected to SDS PAGE followed by western blot with HRP-66005(6*His, His-Tag Antibody) at dilution of 1:10000

Clonality