

# ALOX5 antibody



Catalog Number: 107113

## Product name

ALOX5 antibody

## Specificity

Human, Pig; other species not tested.

## Antibody description

ALOX5 Mouse Monoclonal antibody. Positive WB detected in pig spleen tissue. Observed molecular weight by Western-blot: 78 kDa

## Preparation

This antibody was obtained by immunization of ALOX5 recombinant protein (Accession Number: NM\_000698). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



pig spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:107573(ALOX5 Antibody) at dilution of 1:2000