

CD23 antibody

Catalog Number: 107125

Product name

CD23 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CD23 Mouse Monoclonal antibody. Positive WB detected in Raji cells. Positive IHC detected in human tonsillitis tissue, human lymphoma tissue. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of CD23 recombinant protein (Accession Number: BC064417). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:50-1:500


Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:107125(CD23,FCER2 antibody) at dilution of



1:1000 Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107125(CD23,FCER2 Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded

CD23 antibody



Catalog Number: 107125

human tonsillitis tissue slide using Catalog
No:107125(CD23,FCER2 Antibody) at dilution of

1:200 (under 40x lens). heat mediated antigen
retrieved with Tris-EDTA buffer(pH9).