

CD68 antibody



Catalog Number: 107140

Product name

CD68 antibody

Specificity

Human; other species not tested.

Antibody description

CD68 Mouse Monoclonal antibody. Positive IHC detected in human liver tissue.

Preparation

This antibody was obtained by immunization of CD68 recombinant protein (Accession Number: NM_001251). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human liver tissue slide using (CD68 Antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemistry of paraffin-embedded human liver tissue slide using (CD68 Antibody) at dilution of 1:200 (under 40x lens)