

# CHOP; GADD153 antibody



Catalog Number: 107155

## Product name

Monoclonal

CHOP; GADD153 antibody

## Ig Type

## Specificity

Mouse IgG2a

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

CHOP; GADD153 Mouse Monoclonal antibody.  
Positive WB detected in Tunicamycin treated  
HeLa cells. Observed molecular weight by  
Western-blot: 30 kDa

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of  
CHOP; GADD153 recombinant protein (Accession  
Number: NM\_004083). Purification method:  
Protein A purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Tunicamycin treated HeLa cells were subjected to  
SDS PAGE followed by western blot with Catalog  
No:107155(CHOP; GADD153 Antibody) at dilution  
of 1:1000