

CNPY2, MSAP antibody



Catalog Number: 107167

Product name

CNPY2, MSAP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CNPY2, MSAP Mouse Monoclonal antibody.
Positive IF detected in HepG2 cells. Positive WB detected in HepG2 cells, human placenta tissue, rat liver tissue. Observed molecular weight by Western-blot: 21 kDa

Preparation

This antibody was obtained by immunization of CNPY2, MSAP recombinant protein (Accession Number: NM_014255). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107167(CNPY2,MSAP Antibody) at dilution of 1:1000



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:107167(CNPY2,MSAP Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-mouse IgG(H+L)