

ESR2 antibody



Catalog Number: 107195

Product name

ESR2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ESR2 Mouse Monoclonal antibody. Positive WB detected in MCF-7 cells, human brain tissue, SW 1990 cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of ESR2 recombinant protein (Accession Number: XM_047431078). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107195(ESR2 Antibody) at dilution of 1:2000