

FGF12 antibody



Catalog Number: 107207

Product name

FGF12 antibody

Specificity

Human; other species not tested.

Antibody description

FGF12 Mouse Monoclonal antibody. Positive WB detected in Human brain, human brain tissue. Observed molecular weight by Western-blot: 20 kDa

Preparation

This antibody was obtained by immunization of FGF12 recombinant protein (Accession Number: NM_004113). Purification method: Caprylic Acid/Ammonium Sulfate Precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107207(FGF12 antibody) at dilution of 1:1000