

Fibrinogen beta chain antibody



Catalog Number: 107211

Product name

Fibrinogen beta chain antibody

Specificity

Human; other species not tested.

Antibody description

Fibrinogen beta chain Mouse Monoclonal antibody. Positive WB detected in human heart tissue, fetal human brain tissue. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of Fibrinogen beta chain recombinant protein (Accession Number: NM_005141). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:107211(FGB Antibody) at dilution of 1:1000