

Fibromodulin antibody



Catalog Number: 107213

Product name

Fibromodulin antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

Fibromodulin Mouse Monoclonal antibody. Positive WB detected in Human kidney, human kidney tissue, human stomach tissue. Positive IHC detected in human hepatocirrhosis tissue. Observed molecular weight by Western-blot: 59 kDa

Preparation

This antibody was obtained by immunization of Fibromodulin recombinant protein (Accession Number: NM_002023). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107213(FMOD antibody) at dilution of 1:250



Immunohistochemistry of paraffin-embedded human hepatocirrhosis slide using Catalog No:107213(FMOD Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human hepatocirrhosis slide using Catalog No:107213(FMOD Antibody) at dilution of 1:50