

Flag tag antibody

Catalog Number: 107219

Product name

Flag tag antibody

Specificity

Tag

Antibody description

Flag tag Mouse Monoclonal antibody. Positive IP detected in Recombinant protein protein. Positive WB detected in Recombinant Protein.

Preparation

This antibody was obtained by immunization of Peptide. Purification method: .

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:2000-1:20000

IP: 1:1000-1:10000

Validations



WB result of anti-Flag (Catalog No:117314, 1:8000) with Flag-recombinant protein.



IP Result of anti-Flag tag (IP:Catalog No:117314, 3ug; Detection:Catalog No:117313 1:2000) with Recombinant protein protein lysate 1280ug.