

# Follistatin antibody



Catalog Number: 107222

## Product name

Follistatin antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Follistatin Mouse Monoclonal antibody. Positive IF detected in Hela cells. Positive WB detected in Hepg2 cells, HeLa cells, Raji cells. Observed molecular weight by Western-blot: 32 / 35/ 40 kDa

## Preparation

This antibody was obtained by immunization of Follistatin recombinant protein (Accession Number: NM\_013409). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications

ELISA, WB, IF

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107222(FST antibody) at dilution of 1:1000



Immunofluorescent analysis of Hela cells, using FST antibody Catalog No: at 1:25 dilution and Rhodamine-labeled goat anti-mouse IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).