

FOLR1 antibody



Catalog Number: 107223

Product name

FOLR1 antibody

Immunogen

[Human FOLR1 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

FOLR1 Mouse Monoclonal antibody. Positive WB detected in HeLa cells, MCF-7 cells. Observed molecular weight by Western-blot: 36-38 kDa

Preparation

This antibody was obtained by immunization of FOLR1 recombinant protein (Accession Number: NM_016725). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107223(FOLR1 Antibody) at dilution of 1:1000