

CYP2D6 antibody



Catalog Number: 107226

Product name

CYP2D6 antibody

Immunogen

[Human CYP2D6 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

CYP2D6 Mouse Monoclonal antibody. Positive WB detected in HepG2 cells, fetal human brain tissue. Observed molecular weight by Western-blot: 56 kDa

Preparation

This antibody was obtained by immunization of CYP2D6 recombinant protein (Accession Number: NM_000106). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107226(CYP2D6 Antibody) at dilution of 1:1000