

DAO antibody



Catalog Number: 107235

Product name

DAO antibody

Immunogen

[Human DAO Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

DAO Mouse Monoclonal antibody. Positive WB detected in fetal human brain tissue. Observed molecular weight by Western-blot: 39 kDa

Preparation

This antibody was obtained by immunization of DAO recombinant protein (Accession Number: NM_001917). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107235(DAO Antibody) at dilution of 1:1000