

EEF1D antibody

Catalog Number: 107250

Product name

EEF1D antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

EEF1D Mouse Monoclonal antibody. Positive WB detected in transfection MCF-7 cells. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 31kD+72kD

Preparation

This antibody was obtained by immunization of EEF1D recombinant protein (Accession Number: XM_047421433). Purification method: Caprylic Acid/Ammonium Sulfate Precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200

Validations

Immunohistochemistry of paraffin-embedded human pancreas cancer slide using Catalog No:107250(EEF1D Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human pancreas cancer slide using Catalog No:107250(EEF1D Antibody) at dilution of 1:50