

EIF3D antibody



Catalog Number: 107251

Product name

EIF3D antibody

IHC: 1:20-1:200

IF: 1:50-1:500

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EIF3D Mouse Monoclonal antibody. Positive WB detected in HepG2 cells, A431 cells, HeLa cells, human brain tissue, mouse liver tissue, rat liver tissue. Positive IF detected in HeLa cells. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 66 kDa

Preparation

This antibody was obtained by immunization of EIF3D recombinant protein (Accession Number: NM_003753). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107251(EIF3D antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human brain using Catalog No:107251(EIF3D antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human brain using Catalog No:107251(EIF3D

EIF3D antibody

Catalog Number: 107251

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog No:107251(EIF3D Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)