

EPHA1-specific antibody



Catalog Number: 107255

Product name

EPHA1-specific antibody

Specificity

Human; other species not tested.

Antibody description

EPHA1-specific Mouse Monoclonal antibody. Positive WB detected in HeLa cells, COLO 320 cells, HepG2 cells. Observed molecular weight by Western-blot: 108 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005232). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107255(EPHA1-special Antibody) at dilution of 1:3000