

# HDAC1 antibody



Catalog Number: 107257

## Product name

HDAC1 antibody

## Immunogen

[Human HDAC1 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

HDAC1 Mouse Monoclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 65kD

## Preparation

This antibody was obtained by immunization of HDAC1 recombinant protein (Accession Number: NM\_004964). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:1000-1:10000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107257(HDAC1 antibody) at dilution of 1:4000