

GAMT antibody



Catalog Number: 107288

Product name

GAMT antibody

Specificity

Human, Rat, Mouse, Pig; other species not tested.

Antibody description

GAMT Mouse Monoclonal antibody. Positive WB detected in fetal human brain tissue, mouse liver tissue, pig liver tissue, rat liver tissue. Observed molecular weight by Western-blot: 26 kDa

Preparation

This antibody was obtained by immunization of GAMT recombinant protein (Accession Number: NM_000156). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107288(GAMT Antibody) at dilution of 1:2000