

HRP-conjugated GAPDH antibody



Catalog Number: 107289

Product name

HRP-conjugated GAPDH antibody

Specificity

Human, Mouse, Rat, Zebra Fish; other species not tested.

Antibody description

HRP-conjugated GAPDH Mouse Monoclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 36 kDa

Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM_001289745). Purification method: Protein A Purified.

Formulation

PBS with 0.02% Thimerosal, 50% Glycerol, pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

Validations



Western blot analysis of GAPDH in various tissues and cell lines using Proteintech antibody HRP-60004 at a dilution of 1:10000.