

GEFT antibody

Catalog Number: 107293

Product name

GEFT antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GEFT Mouse Monoclonal antibody. Positive WB detected in mouse skeletal muscle tissue, human brain tissue, mouse brain tissue, rat brain tissue, rat skeletal muscle tissue. Positive IHC detected in human heart tissue, human brain tissue. Observed molecular weight by Western-blot: 53 kDa, 63 kDa

Preparation

This antibody was obtained by immunization of GEFT recombinant protein (Accession Number: NM_133483). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:107293(GEFT antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human heart using Catalog No:107293(GEFT antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:107293(GEFT antibody) at dilution of 1:50 (under 40x lens)