

# GKN1 antibody



Catalog Number: 107298

---

## Product name

GKN1 antibody

## Specificity

human; other species not tested.

## Antibody description

GKN1 Mouse Monoclonal antibody. Positive WB detected in Human stomach. Observed molecular weight by Western-blot: 18-20kd

## Preparation

This antibody was obtained by immunization of GKN1 recombinant protein (Accession Number: BC059778). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



human stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:107298(GKN1 antibody) at dilution of 1:400