

# GLS antibody



Catalog Number: 107299

## Product name

GLS antibody

## Specificity

Human; other species not tested.

## Antibody description

GLS Mouse Monoclonal antibody. Positive WB detected in fetal human brain tissue, HEK-293 cells, rat liver tissue. Observed molecular weight by Western-blot: 58 kDa, 65 kDa

## Preparation

This antibody was obtained by immunization of GLS recombinant protein (Accession Number: NM\_014905). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

## Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107299(GLS Antibody) at dilution of 1:10000