

Glutamine synthetase antibody



Catalog Number: 107301

Product name

Glutamine synthetase antibody

Specificity

Human, Pig; other species not tested.

Antibody description

Glutamine synthetase Mouse Monoclonal antibody. Positive WB detected in Jurkat cells, pig brain tissue, pig liver tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of Glutamine synthetase recombinant protein (Accession Number: NM_002065). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107301(Glutamine synthetase Antibody) at dilution of 1:4000