

# GP73/GOLPH2 antibody



Catalog Number: 107305

## Product name

GP73/GOLPH2 antibody

## Specificity

Human; other species not tested.

## Antibody description

GP73/GOLPH2 Mouse Monoclonal antibody.  
Positive IF detected in HeLa cells. Positive WB  
detected in HEK-293 cells. Observed molecular  
weight by Western-blot: 73 kDa

## Preparation

This antibody was obtained by immunization of  
GP73/GOLPH2 recombinant protein (Accession  
Number: NM\_016548). Purification method:  
Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE  
followed by western blot with Catalog  
No:107305(GOLM1 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol)  
fixed HeLa cells using Catalog No:107305(GOLM1  
Antibody) at dilution of 1:100 and Alexa Fluor 488-  
conjugated AffiniPure Goat Anti-Mouse IgG(H+L)