

GP73/GOLPH2 antibody



Catalog Number: 107305

Product name

GP73/GOLPH2 antibody

Specificity

Human; other species not tested.

Antibody description

GP73/GOLPH2 Mouse Monoclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 73 kDa

Preparation

This antibody was obtained by immunization of GP73/GOLPH2 recombinant protein (Accession Number: NM_016548). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107305(GOLM1 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:107305(GOLM1 Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)